


2519 WOMEN'S HEATHERED POLO

WOMEN'S SIZES XS-3XL

- :: Wicking & Anti-microbial Technology
- :: Peached fabric finish for extra softness
- :: Coordinates with Men's style 3519



Self fabric collar with self-binding for comfort

Open placket

Capped sleeves

Moisture wicking with TopShield™ moisture management

Anti-microbial treatment inhibits the growth of odor-causing bacteria

Fitted styling

Double needle stitches seams for durability

Side-seam vents

2519 WOMEN'S HEATHERED POLO

COMFORT & CONVENIENCE

- :: Made exclusively with TopShield™ moisture management treatment. This hi-tech treatment applied to the fabric wicks moisture away from your body, keeping you comfortable & dry
- :: Antimicrobial treatment inhibits the growth of odor-causing bacteria
- :: 100% heathered polyester jersey with a peached finish for extra softness
- :: Tagless neck labeling and self-binding at collar prevents chafing against skin
- :: Self fabric collar with self-binding for comfort
- :: Open placket, capped sleeves, side-seam vents and fitted styling
- :: Suitable for corporate sponsored events, school & team wear, everyday wear

FABRIC

100% heathered polyester jersey 3.83 oz/yd2

SIZES

Women's sizes XS-3XL

COLORS

057 Orchid Heather, 208 Graphite Heather, 240 Navy Heather, 176 Red Heather



057 ORCHID HEATHER



208 GRAPHITE HEATHER



240 NAVY HEATHER



176 RED HEATHER

HEATHERED POLO WOMEN'S 2519	XS	S	M	L	XL	2XL	3XL	TOLERANCE
BODY LENGTH measure from neck (cb neck)	24½	25	25½	26¼	27¼	28	29	¾ +/-
ACROSS SHOULDER measure from back	14	14½	15	15¾	16⅝	17¾	18⅝	¼ +/-
CHEST WIDTH 1 inch below armhole	17¼	18¼	19¼	21	23	25¼	28	¾ +/-
WAIST WIDTH	16½	17½	18½	20¼	22½	25	27½	¾ +/-
BOTTOM OPENING WIDTH RELAXED 1/2	18½	19½	20½	22¼	24¼	26¾	29¼	¾ +/-
SLEEVE LENGTH FROM CB short	13	13½	14	15	16¼	17	17¾	¼ +/-

KEY TO COMMONLY ABBREVIATED SPEC INFORMATION (ALL MEASUREMENTS ARE IN INCHES)

CB = CENTER BACK

CHEST 1 INCH DOWN = MEASURE THE FRONT OF THE GARMENT 1 INCH DOWN FROM THE ARMHOLE

1/2 = MEASURE ONLY THE FRONT OF THE GARMENT, THEN DOUBLE THAT MEASUREMENT TO GET THE TOTAL CIRCUMFERENCE

RELAXED = DO NOT STRETCH THE GARMENT IF IT HAS ELASTIC TO MEASURE, SIMPLY LEAVE THE ELASTIC RELAXED

STRETCHED = STRETCH THE GARMENT IF IT HAS ELASTIC TO MEASURE, MEASURE STRETCHED