

# 2811 WOMEN'S CLASSIC WICKING POLO

WOMEN'S SIZES XS-3XL

- :: Wicking & Anti-microbial Technology
- :: Provides UV protection
- :: Coordinates with Men's style 3811



Taped neckline for clean, finished appearance

Solid body and rib knit collar accented with white tipping

5-button placket with dyed-to-match pearlized buttons

Set-in sleeve design

Moisture wicking with TopShield™ moisture management

Fitted styling

Anti-microbial treatment inhibits the growth of odor-causing bacteria

100% polyester piqué knit



BACK VIEW

**COMFORT & CONVENIENCE**

- :: Made exclusively with TopShield™ moisture management treatment. This hi-tech treatment applied to the fabric wicks moisture away from your body, keeping you comfortable and dry
- :: Anti-microbial treatment inhibits the growth of odor causing bacteria, leaving you feeling fresh
- :: Provides UV protection (tested and rated to provide a UPF of up to 27).
- :: Tagless neck labeling and self-binding at the collar prevents chafing against skin
- :: Extra button conveniently located on inner side seam label
- :: Suitable for corporate sponsored events, school & team wear, everyday wear
- :: Coordinates with Men's style 3811

**FABRIC**

100% polyester piqué knit 4.87 oz/yd2 (8.11 oz/lyd, 165gsm)

**SIZES**

Women's Sizes XS-3XL

**COLORS**

010 Black, 040 Navy, 060 Red, 070 Royal, 080 White, 180 Pink

**GARMENT CARE**

Machine wash cold separately. Do not use fabric softeners. Only non-chlorine bleach when needed. Tumble dry low. Cool iron if needed.



CLASSIC WICKING POLO WOMEN'S 2811	XS	S	M	L	XL	2XL	3XL	TOLERANCE
<b>BODY LENGTH</b> measure from neck (cb neck)	23¾	24¼	24¾	25½	26¼	27	27¾	¾ +/-
<b>ACROSS SHOULDER</b> imaginary shoulder points	14½	15	15½	16¼	17	18	19	¼ +/-
<b>BODY WIDTH</b>	16¾	17¾	18¾	20¼	21¾	23¾	25¾	¾ +/-
<b>BOTTOM OPENING WIDTH RELAXED 1/2</b>	18	19	20	21½	23	25	27	¾ +/-

KEY TO COMMONLY ABBREVIATED SPEC INFORMATION (ALL MEASUREMENTS ARE IN INCHES)

CB = CENTER BACK

CHEST 1 INCH DOWN = MEASURE THE FRONT OF THE GARMENT 1 INCH DOWN FROM THE ARMHOLE

1/2 = MEASURE ONLY THE FRONT OF THE GARMENT, THEN DOUBLE THAT MEASUREMENT TO GET THE TOTAL CIRCUMFERENCE

RELAXED = DO NOT STRETCH THE GARMENT IF IT HAS ELASTIC TO MEASURE, SIMPLY LEAVE THE ELASTIC RELAXED

STRETCHED = STRETCH THE GARMENT IF IT HAS ELASTIC TO MEASURE, MEASURE STRETCHED