# **5535 WOMEN'S RADIUS QUILTED VEST**

## **WOMEN'S SIZES XS-3XL**

- :: 100% dobby polyester shell made exclusively from recycled PET bottles
- :: Coordinates with men's style 9535







### **COMFORT & CONVENIENCE**

- :: Ultra-lightweight & packable
- :: 100% dobby polyester shell made exclusively from recycled PET bottles
- :: Lined with PrimaLoft® SPORT for premium insulation performance
- :: Engineered for superior warmth
- :: Ultra compressible fabric allows vest to pack into side pocket
- :: Features narrow, horizontal quilted pattern with contrast-colored zippers & lining
- :: Full zip styling with wind flap
- :: Zippered on-seam pockets
- :: Dropped tail hem for added coverage
- :: Suitable for hospital nurses, employee gifts, corporate wear, or layering when on-the-go
- :: Coordinates with Men's Radius Quilted Vest, style 9535

**FABRIC** 

100% dobby polyester lined with PrimaLoft® SPORT ()60 gm²)

**SIZES** 

XS-3XL

#### **COLORS**

006 Black/Grey, 096 Navy/Grey, 197 Cobalt/Grey, 296 Grey/Pink

#### **GARMENT CARE**

Machine wash separately cold water inside out. Gentle cycle with like colors. Do not bleach. Do not use fabric softeners. Line dry. Do not iron.

006 BLACK/GREY

096 NAVY/GREY

197 COBALT/GREY

296 GREY/PINK



RADIUS QUILTED VEST WOMEN'S 5535	XS	S	М	L	XL	2XL	3XL	TOLERANCE
BODY LENGTH measure from CB neck	241/2	25	25½	261/4	27	275/8	281/4	1/2 +/-
ACROSS SHOULDER measure from back	141/4	14¾	151/4	16	16%	18	191/8	1/2 +/-
CHEST WIDTH 1" below armhole (following curve to include chest depth)	18½	19½	20½	221/4	241/4	26¾	291/4	1/2 +/-
WAIST WIDTH	17	18	19	20¾	22¾	251/4	27¾	1/2 +/-
BOTTOM OPENING WIDTH relaxed	191/4	201/4	211/4	23	25	271/2	30	1/2 +/-

KEY TO COMMONLY ABBREVIATED SPEC INFORMATION (ALL MEASUREMENTS ARE IN INCHES)

CB = CENTER BACK

CHEST 1 INCH DOWN = MEASURE THE FRONT OF THE GARMENT 1 INCH DOWN FROM THE ARMHOLE

1/2 = MEASURE ONLY THE FRONT OF THE GARMENT, THEN DOUBLE THAT MEASUREMENT TO GET THE TOTAL CIRCUMFERENCE
RELAXED = DO NOT STRETCH THE GARMENT IF IT HAS ELASTIC TO MEASURE, SIMPLY LEAVE THE ELASTIC RELAXED

STRETCHED = STRETCH THE GARMENT IF IT HAS ELASTIC TO MEASURE, MEASURE STRETCHED

